



Tunstall

Introducing

Tunstall GO™

Nick Brown
Business Development Director
South East Region

Challenges of post pandemic anxiety

- “Figures show that the equivalent of 19 million adults in Great Britain report high levels of anxiety.”
- "Prior to the coronavirus pandemic we consistently saw lower anxiety ratings in those aged 65 years and older, but now we are seeing the highest levels of anxiety amongst this group".
- *Dawn Snape, Assistant Director of Sustainability and Inequalities Division, Office for National Statistics*

How do we support people post lockdown?

- 18 year old wheelchair user
 - Previously independent in the community
 - Panic attacks
 - Fear of leaving the house
- 83 year old
 - No major health issues
 - Reluctant to go out in the community
 - Palpitations



Overview



LOCATION



SUPERB SPEAKERPHONE



AUDIBLE NOTIFICATIONS



REMOTE UPDATES

FEATURES

- Fall detection
- 1-2 month battery
- Wireless charger
- Illuminating notification ring around device
- System check function for battery, coverage and GPS
- Flight mode
- Installer mode
- IP67
- Hypoallergenic chain



SmartCare Locate

The screenshot shows the SmartCare Locate interface. On the left is a table of events, and on the right is a map of a coastal area with a red pin indicating a location. A pop-up window provides details for the selected event.


Account (Code) (Method)	Date Time
5980 (Lost) (GPS Rough)	10/May/2018, 10:48:37 am
20000154 (Check) (GPS Rough)	10/May/2018, 9:45:50 am
20000155 (Lost) (GPS Rough)	10/May/2018, 4:35:47 am
5984 (Check) (GPS Rough)	10/May/2018, 2:52:18 am
5985 (Help) (GPS)	10/May/2018, 1:24:16 am
-5985 (1100) (GPS)	10/May/2018, 1:24:16 am
-5985 (1100) (GPS)	10/May/2018, 1:19:44 am
-5985 (1100) (GPS)	10/May/2018, 1:14:15 am
-5985 (1100) (GPS)	10/May/2018, 1:09:14 am
-5985 (1100) (GPS)	10/May/2018, 1:05:14 am
-5985 (1100) (GPS)	09/May/2018, 11:35:41 pm
-5985 (1601) (GPS)	09/May/2018, 11:32:40 pm
20000152 (Check) (GPS Rough)	10/May/2018, 1:21:06 am
5972 (Check) (GPS Rough)	09/May/2018, 9:00:02 pm
20000151 (Check) (GPS Rough)	09/May/2018, 8:50:46 pm

Map Pop-up Details:

- Account: 5985
- EventCode: 1100 (Help Activation)
- Serial Number: 20000204
- Method: G (GPS)
- Latitude: -43.558861875
- Longitude: 172.70662375
- Cellular Coverage: -87dBm (Excellent)
- Date Time Reported: 5/10/2018 1:24:16 AM
- Data: 16.0

- GPS coordinates are transmitted via cellular data to SmartCare® Locate
- Direct SmartCare Locate URL sent with alarm transmission

Battery and charging



The GO is provided with an inductive charger for ease of use



Rechargeable lithium-ion battery approximately 1 -2 months before requiring a charge



Voice messages and lights indicate at approximately 20% capacity and is reported to monitoring



Once charging to an acceptable level is completed a report will be sent to monitoring



