



# Vibe

Sensory and Communication Aid



Creating products that speak for themselves...

[we-cair.com](http://we-cair.com)



# Vibe


Sensory and Communication Aid

**The Vibe is a compact, all-in-one sensory and communication aid** that promotes reassurance, independence, and peace of mind for people with diverse support needs. It responds to connected sensors, such as smoke alarms and doorbells, using a fully customisable mix of visual, audible, vibratory, and spoken alerts.

With bright strobe lights, LED icons, clear tones, and voice prompts, every alert is unmissable. A vibrating pillow pad offers discreet night-time notifications, and voice-guided setup makes installation quick and tool-free.

Five LED icons (including dedicated fire and gas indicators) and programmable locations help users identify each alert. Mains powered with a 24-hour backup battery, the Vibe supports a wide range of sensors, including movement, medication, bed occupancy, and door/window contacts.

Ideal for people with sensory needs, autism, dementia, or learning disabilities, the Vibe can also notify carers when telecare sensors are triggered.



Creating products that speak for themselves...

[we-cair.com](http://we-cair.com)



# Vibe

Sensory and Communication Aid

**There are endless possibilities with the Vibe, which can be used in a multitude of ways with a huge range of wireless peripherals.**

From visually alerting someone living with hearing loss that a smoke detector has activated, or alerting someone with visual needs that there is a flood in the kitchen, or waking a carer by the vibrating pillow alert that someone needs support at night; the Vibe is the solution.



Creating products that speak for themselves...

[we-cair.com](http://we-cair.com)



# Vibe

Sensory and Communication Aid

Multiple Wireless  
Sensor Inputs

520Hz Low  
Frequency Tone

24Hr Battery Backup


2 Pre-Programmed &  
3 Programmable LEDs

Clear Speech Alerts

Bright Strobe Light

Outputs for 2 Wired  
Vibrating Pillow Pads

Clear Hi Definition Audio



Creating products that speak for themselves...

[we-cair.com](http://we-cair.com)



# Vibe

Sensory and Communication Aid



Creating products that speak for themselves...

[we-cair.com](http://we-cair.com)



# Vibe

Sensory and Communication Aid

The Vibe can alert a user when it's time for them to take their medication. As well as a dedicated voice announcement, visual and vibratory alerts can ensure medication from the pill dispenser is never missed again..

*It's time to take your medication!*

**Vibe**  
Sensory & Communication Aid

**Pill+**  
Medication Dispenser



Creating products that speak for themselves...

Medicine Prompts using Vibe and Pill +





# Vibe

Sensory and Communication Aid

When an environmental sensor, such as a smoke alarm is activated, the user can be woken up by the vibrating pillow pads that are connected to the Vibe, giving them additional peace of mind..

*FIRE, FIRE, FIRE!*

Vibe  
Sensory & Communication Aid



Creating products that speak for themselves...

Environmental Alerts using Vibe and Smoke, Heat & CO





# Vibe

Sensory and Communication Aid

**Flood**  
Flood Sensor

**Vibe**  
Sensory & Communication Aid

*Flood Detected!*

When a potential flood goes unnoticed, the damage and risk can be significant. Ensure users are alerted immediately if a flood is detected, so that action can be taken straight away.



Creating products that speak for themselves...

Flood Alerts using Vibe and Flood





# Vibe

Sensory and Communication Aid

Contact  
Door/Window Sensor

Door or Window Open!



Vibe  
Sensory & Communication Aid

Whether it's for a user who wants to be alerted of visitors coming into their home, or a carer with sensory loss who needs to be alerted when the person they're supporting has opened a door or window, the Vibe can do it all.



Creating products that speak for themselves...

Door/Window Open/Closed Alerts using Vibe and Contact






# Vibe


Sensory and Communication Aid

**Move**  
Motion Sensor



The **Move** is designed to detect motion. When used with the Vibe, it can notify a carer when movement has been detected. It can also offer additional security for those living alone if movement is detected in their home.

**Vibe**  
Sensory & Communication Aid



**Motion Detected!**



Creating products that speak for themselves...

Motion Alerts using Vibe and Move





# Vibe

Sensory and Communication Aid

The **Chime** can replace a traditional doorbell, and be used to notify the user when someone activates it. Adapt the Vibe to notify the user in a way that works for them, ensuring they never miss a visitor again.



There is someone at your door!

Vibe  
Sensory & Communication Aid

Chime  
Doorbell/Bogus Caller/Place Anywhere Button



Creating products that speak for themselves...

Push Button Alerts using Vibe and Chime





# Vibe

Sensory and Communication Aid

The **Request** can be used as a panic button and can be activated by a person or support animal to alert a carer that support is needed. Once activated, the Vibe will immediately notify the user.

**Request**  
Large Button Radio Alarm


*Panic Button Activated!*

**Vibe**  
Sensory & Communication Aid



Creating products that speak for themselves...


Push Button Alerts using Vibe and Request





# Vibe

Sensory and Communication Aid



Creating products that speak for themselves...

[we-cair.com](http://we-cair.com)



# Vibe

Sensory and Communication Aid



# THANK YOU

Lowri Jones

[lowri@we-cair.com](mailto:lowri@we-cair.com)



Creating products that speak for themselves...

[we-cair.com](http://we-cair.com)